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## BAS2 / FOR2 / LIT2

- Aufbereitungsanleitung
- Processing instructions
- Instructions de retraitement
- Instrucciones de acondicionamiento
- Istruzioni di trattamento
- Instructie voor het opnieuw gebruiksklaar maken
- Instruções de reprocessamento
- Rengöringshandledning
- Oparbejdningsvejledning
- Οδηγίες επανεπεξεργασίας
- Hazırlık talimatları
- Uputa za pripravu
- Návod na přípravu
- Taastötlemise juhend
- Előkészítési útmutató
- Sagatavošanas instrukcija
- Apdrojimo instrukcija
- Instrukcja reprocesowania
- Navodila za pripravo za uporabo

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## Specifications for the reprocessing of endoscopic instruments

In order to guarantee functional reliability and a long service life, reusable endoscopic instruments must be processed after every examination in accordance with the applicable guidelines and standards. These specifications are based on DIN EN ISO 17664 and "Hygiene requirements for the reprocessing of flexible endoscopes and endoscopic additional instruments" published by the Robert Koch Institute (RKI).

The automatic reprocessing process in the washer-disinfector (WD) for reusable medwork products as described below has been validated by medwork and constitutes the specifications for the operator. Please note that the operator is solely responsible for carrying out and monitoring the correct and proper reprocessing.

## Products

The specifications published here contain important information for the reliable and effective reprocessing of medwork products and apply in conjunction with the respective instructions for use of the products.

The type and scope of reprocessing depends on the application of the medical device. Therefore, the operator is responsible for the correct classification of the medical devices and as a result for specifying the type and scope of reprocessing (see KRINKO/BfArM recommendation, Section 1.2.1 Risk assessment and classification of medical devices before reprocessing). It must particularly be ensured that standardised and validated processes are applied. If reusable medwork products are used for patients suffering from a variant of Creutzfeld-Jakob disease, they must not be processed after use. The particular recommendations of the RKI are applicable here.

medwork products that are labelled as single use devices (⊗-symbol) may not be prepared or resterilised for reuse. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health, or result in illness, injury or death. Reuse, reprocessing and resterilisation also present the additional risk of patient or product contamination, as well as the risk of cross contamination including the transmission of infectious diseases. Product contamination can result in patient illness, injury or death. All instruments must be cleaned, disinfected and sterilised before every use. This is particularly important when using an instrument for the first time after delivery, because all instruments are supplied in sterile condition (cleaning and disinfection after removal of protective shipping packaging;

sterilisation after packing). Particular care is required for cleaning instruments with hollow bodies. A special rack with provision for rinsing cavities is required for machine cleaning and disinfection of such instruments, e.g. racks for micro-invasive instruments (MIC trolley). The equipment of these racks varies depending on the manufacturer. The cleaning performance depends on the equipment and adaptability of the instruments to the rack and must be confirmed by validation of the processes. Please follow the manufacturer's information on the compatibility of materials for cleaning, disinfection and sterilisation. For ultrasonic cleaning use only cleaners approved for the process and maintain the temperature range recommended by the manufacturer.

## Premises

The reprocessing must only be conducted in suitable rooms specially equipped for this purpose and that comply with the requirements of the Committee for Hospital Hygiene and Infection Prevention of the RKI.

## Personnel

Please note that additional endoscopic instruments must only be processed by trained persons (Category IV) who have the necessary training, level of knowledge and experience for this work. You should wear a face mask, safety glasses, an operating gown impervious to liquids and cut-resistant gloves for protection from potential contamination with pathogenic germs and potentially hazardous cleaning chemicals.

## Machine reprocessing

### Equipment for precleaning and manual cleaning:

Cleaning agent: non-protein-fixing VAH-listed instrument disinfectant with cleaning action, e.g. based on quaternary compounds,

guanidine/guanidine derivatives, gigasept AF forte® (Schülke & Mayr GmbH) is validated

Plastic brushes	1st brush: brush head diameter 2.5 mm 2nd brush: brush head length 100 mm, brush head diameter: 7.0 mm
Disposable cloth or sponge	lint-free disposable cloth or sponge, validated: B. Braun, article number: 19726
Disposable syringe:	20 ml syringe, validated: B. Braun, article number: 4606205V
Ultrasonic bath:	validated: Bandelin, Sonorex RK1028
Mains water:	running water, temperature 20 +/- 2°C, minimum potable water quality
Sink:	sink for mains water
Transport box:	lockable transport container

### Step 1: Precleaning:

Clean used instruments as soon as possible after endoscopy with a lint-free disposable cloth or sponge in order to minimise drying on of organic material and chemical residues as much as possible and prevent contamination of the environment. For the protection of personnel and the environment, transport the contaminated instrument only in a lockable container from the examination room to the conditioning room. We recommend reprocessing of the instruments as soon as possible, but in any case no later than 2 hours after use. Temporary storage of used instruments with impurities such as blood residues may cause corrosion.

### Step 2: Manual cleaning:

#### ATTENTION! Special Instructions

The instruments must be disassembled for cleaning and disinfection depending on the device. Follow the instructions for disassembly in the instructions for use. Do not use any solutions containing aldehydes as these can fix proteins and thus adversely affect the success of the cleaning.

Prepare the cleaning solution (e.g.: gigasept AF forte®2%) as directed by the manufacturer and place it in an ultrasonic bath. If possible, place the instruments immediately after use in the cleaning solution prepared strictly in accordance with the concentrations, interactions and exposure times specified by the manufacturer. Carry out all following cleaning steps below the surface of the liquid in order to prevent spraying of contaminated liquid. Use only cleaning solutions whose suitability and effectiveness have been documented and declared accordingly in the expert reports of the reprocessing and equipment manufacturers. The cleaning solution must be replaced at least daily and immediately in the case of visible contamination. Insert the instrument with a diameter of 20 cm into the basket. Wipe the wires thoroughly with a lint-free cloth soaked in cleaning solution. Finally, brush the wires for at least 1 minute (1 minute has been validated) with the cleaning brush (validated brush head length: 100 mm, brush head diameter: 7.0 mm). During the soaking time brush the parts of the immersed instrument that are difficult to access with a soft brush (validated brush head diameter 2.5 mm). Pay particular attention to the critical, difficult to access parts where visual inspection of the cleaning effect is not possible. Brush the Luer-Lock attachment and the cavities with a soft brush (validated brush head diameter 2.5 mm). Insert the brush head into the lumen, push the brush carefully through the lumen and then pull it out carefully. Finally, push the brush into the lumen from the distal end and pull it out carefully.

Then rinse the lumen once with the cleaning solution using a 20 ml syringe. Fill the lumen with cleaning solution with no air bubbles. Run the ultrasonic bath with the instruments for at least 5 minutes (5 minutes has been validated). Then rinse the lumen at least once with cleaning solution using a 20 ml syringe to remove dissolved impurities. The total time in the cleaning solution must be at least 15 minutes. After that remove the instrument from the cleaning solution and place it in a sink with water (minimum potable water quality) for at least 1 minute. Rinse the lumen at least once with cleaning solution using a 20 ml syringe to flush out completely the cleaning solution residue.

## ATTENTION! Special Instructions

Please inspect the instrument after manual conditioning to make sure that it is clean. If you can still see visible impurities, please repeat the above procedures.

## Equipment for machine cleaning and disinfection in the WD (washer-disinfector):

Washer-disinfector:	washer-disinfector in accordance with DIN EN ISO 15883-1/-2 with thermal programme (temperature 90°C to 95°C)
Cleaner:	mild alkaline cleaner, validated: Ecolab, Sekumatic® MultiClean

## Step 3: Machine cleaning and disinfection

Lay or position the individual parts in a suitable mesh tray or rack so all internal and external surfaces of the instruments can be cleaned and disinfected. Insert instruments with a diameter of 20 cm. The instruments must be attached to the rack (e.g. MIC device trolley) to ensure that the internal lumens are rinsed. The attachment must be checked to ensure that it is tight before starting the machine and at the end of the process. Close the washer-disinfector and start the programme.

See the table below for the programme:

Programme step	Water	Dosage	Time	Temperature
Prewash	CW		5 min	
Cleaner dosage		As per manufacturer's directions		As per manufacturer's directions
Clean	DI		10 min	55°C
Rinse	DI		2 min	
Disinfect	DI		3 min	A <sub>v</sub> value > 3000- e.g. 90°C, 5 min)
Dry			15 min	Up to 120°C

1 CW = cold water

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3 Regulators may specify different processing regulations in their areas of responsibility (parameters for disinfection performance).

Remove all medical devices at the end of the programme. Check that the load is dry and if necessary dry the load with medical compressed air. After emptying the WD inspect the load to make sure that it is clean. If contamination is still visible, steps 2 and 3 of the reprocessing instructions must be repeated.

## Equipment for maintenance, inspection and testing

Oil spray: silicone-free, steam-sterilisable oil spray approved for use with medical devices, e.g. Sterilit®, oil spray supplied by B. Braun Aesculap

## Step 4: Maintenance, inspection and testing

Visually inspect the instrument to check that it is clean, undamaged and functional, using a lighted magnifier if necessary (3.6 dioptre). Test the function of the instrument and inspect it for damage and wear. Pay particular attention here to functional movement and completeness of all components. Where applicable, use suitable oil spray on pivot points of moving parts.

## ATTENTION! Special Instructions

Do not use any products which show kinks, cracks, fractures, discoloration, surface changes or similar.

## Equipment for packaging

Foil-paper packaging: validated: steriCLIN, art. no. 3FKFS23011 and 3FKFS230114

Sealing unit: validated: HAWO, type 880 DC-V

## Step 5: Packing

Push the traction cable into the tube until the distal end of the instrument still projects from the tube. Do not coil the instrument any tighter than a diameter of 20 cm and place it in the packaging together with the small parts. A suitable process (sterile barrier system in accordance with DIN EN ISO 11607 (single)) must be followed to pack the instrument.

## ATTENTION! Special Instructions

The instruments can be packed in packages of one or two instruments. The packages must be large enough so the sealing seam is not under tension. (In the packaging process described here two instruments are packed in one foil-paper package.) After heat-sealing the sealing seam must be visually checked for any faults. If any faults are found, the package must be opened and the instrument packed and sealed again.

## Equipment and process for sterilisation:

Steriliser: steriliser in accordance with DIN EN 285 or small steam steriliser in accordance with DIN EN 13060, Type B process  
Process: steam sterilisation with fractionated pre-vacuum, 134°C, retention time min. 5 min (longer retentions times are possible), (validated 134°C, 5 minutes).

## Step 6: Sterilisation

Place the packed devices in the sterilising chamber and start the programme. At the end of the programme remove the instruments and allow them to cool. Then inspect the packages for damage and penetration of moisture. Suspect packages must be evaluated as unsterile. The device must be packed and sterilised again.

## Step 7: Storage

Storage location (protected from dust and moisture as a minimum) and storage time with the user according to the specifications.

## Reusability

The instruments can be reused ten times assuming appropriate care and that they are uncontaminated. Any additional reuse or the use of damaged and/or contaminated instruments is within the responsibility of the operator. No liability will be accepted in case of misuse.

## Additional information:

Further information on conditioning can be found in the internet at:

[www.rki.de](http://www.rki.de)

[www.esgema.org](http://www.esgema.org)

## List of all icons used on medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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Čistilo: srednje alkalno čistilo, potrejeno je: podj. Ecolab, Sekumatic © MultiClean