FOR-Series

FOR2-H3

FOR2-J3

CE

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Instructions for use

Intended use

The FOR Series instruments are used for the manipulation and extraction of foreign bodies from the upper and lower gastrointestinal tract.

Product features


Figure 1:

Handle FOR2-H3-23-230
and FOR2-J3-23-230



Figure 2:

Disassembly and assembly

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

All instruments labelled as unsterile must be prepared in accordance with the applicable medwork preparation instructions before use for the first time.

All medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

Indications

Resected tissue (e.g. polyps), food bolus impaction, foreign body ingestion.

Contraindications

The contraindications for the FOR series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy, rectoscopy and enteroscopy.

Possible complications

Perforation, haemorrhage, infection, sepsis, allergic reaction to medicines, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Attention! For the extraction of sharp objects we recommend using an endoscope cap or an overtube. Deploy the foreign body retrieval instrument only when visible at the distal end of the instrument. During foreign body retrieval in the upper gastrointestinal tract make sure the foreign body is completely and firmly anchored in the foreign body retrieval instrument so that any loss of the foreign body will not lead to aspiration.

Instructions for use

Perform the diagnostic endoscopy and adjust the tip of the endoscope to the foreign body to be extracted. Guide the foreign body removal instrument with closed grippers through the working channel of the endoscope until the tip of the instrument comes into endoscopic view. Open the gripper by operation of the handle and grasp the foreign body for extraction. Close the grippers. Carefully withdraw the endoscope with the foreign body. Pull the instrument out of the endoscope slowly to avoid contamination of the patient and users with blood and other body fluids and the associated risk of infection.

After completion of the surgical procedure

Reusable instruments must be correctly prepared. Follow the enclosed instructions for reprocessing.

If the instruments can not be used anymore, they must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

Disassembly and assembly

Open by turning it and remove the clamping grip from the pulling rod for conditioning the FOR2-H3-23-230 and FOR2-J3-23-230. Unscrew the screw cap and remove the internal rubber seal. Completely remove the traction cable with the gripper from the tube. Send all parts for reprocessing.

Take to prevent contamination during disassembly.

After conditioning and before starting use, assemble the FOR2-H3-23-230 or FOR2-J3-23-230 instruments by inserting the traction cable from the proximal end to the distal end of the tube. Place the rubber seal and the screw cap on the traction cable (for the sequence see Figure 2). Then attach the clamping grip by screwing it onto the pulling rod.

List of all icons used on medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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