
BAS-Series _ ROCKSTAR

BAS3-K2

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
Intended use

The BAS series instruments are used for endoscopic removal of stones from the bile duct and pancreatic duct.

Product features

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

ⓧ ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

Choledocholithiasis.

Contraindications

The contraindications for the BAS series correspond to the specific contraindications for endoscopic retrograde cholangiopancreatography.

Possible complications

Perforation, cholangitis, haemorrhage, aspiration, infection, sepsis, allergic reaction to contrast agents or medicines, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, inflammation, pressure necroses, foreign body in patient, injury e.g. bruise, obstruction, occlusion, pneumonia, thrombosis, tissue damage, hematoma, pain.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious

incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

When using the product, observe the guide wire diameter recommended on the product label.

When extracting stones from the duct system, it may be necessary to perform a papillotomy if the stone will not pass the papilla. Observe all safety measures, warning notices and contraindications when doing this.

ATTENTION! For stone extraction baskets of the BAS series, only system components of the LIT series may be used. FUJIFILM medwork™ assumes no product liability if the product is used with other lithotripsy system components.

Instructions for use

Insert the retracted stone extractor into the working channel of the endoscope and push it forwards a little at a time until the tube tip is visible in the endoscopic image. If access to the bile duct is already provided via a guide wire, insert the end of the guide wire into the guide wire sheath and push the instrument forwards until the tip of the instrument is visible in the endoscopic image. Using the guide wire if necessary, probe for the required ductal system until you have reached the area behind the stone, and capture the stone carefully under radiological guidance. The stone extractor is retracted and extended by operation of the dual-ring handle. The contrast agent can be introduced using the Luer-Lock attachment. Small stones can be extracted by withdrawing the tube via the papilla. If the stone cannot pass through the papilla, the instrument can also be used for mechanical lithotripsy. In order to carry out lithotripsy on a stone, the adjustment wheel must be turned against the dual-ring handle in order to achieve the level of force required for lithotripsy. Following lithotripsy, the adjustment wheel must be turned back to its starting position in order to be able to open the basket. Slowly pull the instrument out of the endoscope to prevent contamination of the patient and the user with blood and other bodily fluids, and an associated risk of infection.

Non-endoscopic mechanical lithotripsy under radiological guidance:

Once the handle has been separated (forceps) from the tube, the endoscope is removed. The lithotripter spring is pressed forward along the traction wire until the stone is reached, and connected with the lithotripter handle. In this regard, please also observe the instructions for use provided for the LIT series. Use the X-ray image to ensure that the distal end of the spiral is placed on one axis with the basket before you start the lithotripsy. By turning the tension screw on the handle, the pressure applied to the stone is gradually increased. Pause for a few seconds before continuing to turn, as the tension is often already sufficient to crush the stone. The BAS series stone extractors that support lithotripsy feature a defined, predetermined breaking point. Once crushing has been completed, the basket with the lithotripter spring is carefully removed from the duct. The stone calculus can be retrieved from the duct using a stone extraction balloon or a stone extractor.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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