BIO-Series _ EndoBite²

BI01-C3


BI01-C4

BI01-D3

BI01-D4

BI01-G4

BI01-H4

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
Intended use

The BIO series instruments are used for taking tissue samples from the upper and lower intestinal tract and from the bronchial tract.

Product features**Qualification of user**

The instruments must only be used with comprehensive knowledge of the technical principles of gastro-intestinal or bronchial endoscopes. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

Polyps, unclear lesions of the gastrointestinal and bronchial mucosa.

Contraindications

The contraindications for the BIO series correspond to the specific contraindications for oesophagogastroduodenoscopy, endoscopic retrograde cholangiopancreatography and bronchoscopy.

Possible complications

Perforation, haemorrhage, aspiration, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, pneumothorax, tissue damage.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or

malfuction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Attention! The biopsy forceps must be closed during insertion, advance and retrieval from the endoscope, otherwise the working channel of the endoscope may be damaged. Tissue samples should only be taken under direct endoscopic observation.

Instructions for use

Perform the diagnostic endoscopy and adjust the endoscope tip to the sample area. Advance the biopsy forceps, with the branches closed, through the working channel of the endoscope until the tip of the biopsy forceps comes into endoscopic view. Open the branches by operating the control and apply the open biopsy forceps to the mucous membrane. The central stylet prevents slippage on firm tissue. Close the branches carefully with light pressure by operating the control. High pressure is not necessary for cutting the tissue and may damage the internal traction cable of the biopsy forceps. Pull the biopsy forceps with closed branches out of the working channel of the endoscope and prepare the tissue sample for testing.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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