**FUJIFILM**

**medwork**

## **FOR-Series \_ OCTOPUS**

**FOR1-N1**

**FOR1-N2**

**FOR1-N3**

de en fr es it nl pt sv da el tr  
hr cs et hu lv lt pl sl bg ga mt  
ro sk fi no is

404450305797

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
## Instructions for use

### Intended use

The FOR Series instruments are used for the manipulation and extraction of foreign bodies from the upper and lower gastrointestinal tract.

### Product features

### Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

### General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

**ATTENTION!** Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

### Indications

Resected tissue (e.g. polyps), food bolus impaction, foreign body ingestion.

### Contraindications

The contraindications for the FOR series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy, rectoscopy and enteroscopy.

### Possible complications

Perforation, haemorrhage, infection, sepsis, allergic reaction to medicines, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, foreign body in patient, abrasion, laceration.

### Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field

service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Attention!** For the extraction of sharp objects we recommend using an endoscope cap or an overtube. Deploy the foreign body retrieval net only under vision of the distal end of the instrument. During foreign body retrieval in the upper gastrointestinal tract make sure the foreign body is completely and firmly anchored in the foreign body retrieval net so that any loss of the foreign body will not lead to aspiration.

### Instructions for use

Perform the diagnostic endoscopy and adjust the tip of the endoscope to the foreign body to be extracted. Advance the closed foreign body retrieval net through the working channel of the endoscope until the distal tip of the tube comes into endoscopic view. Guide the distal tip of the tube past the foreign body. Open the net completely and place the net over the foreign body to be retrieved by operation of the dual-ring handle. For optimal placement of the rotatable foreign body retrieval net, turn the handle to facilitate the retrieval of foreign bodies that are difficult to reach. Close the net by operating the dual-ring handle. Pull the foreign body towards the endoscope under endoscopic vision, giving the foreign body a continuous slight pull. Carefully withdraw the endoscope from the body together with the foreign body collected. Following extubation of the endoscope the foreign body can be removed from the foreign body retrieval net.

### After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

### List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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