

medw|ork®
for mindful people

#clean

Rounded tip protects the working channel

Compact bristle structure for maximized cleaning efficiency

Optimal force transmission during insertion

RACCOON Brush



Product Code	Order No.	Characteristics	ø Brush	Length	Unit
⊗ CLE1-B1-30-16-100	500458	For bronchoscope, 1 brush	3,0 mm	100 cm	50
⊗ CLE1-B2-30-16-100	500459	For bronchoscope, 2 brushes	3,0 mm	100 cm	50
⊗ CLE1-G1-30-16-180	500468	For gastroscope, 1 brush	3,0 mm	180 cm	50
⊗ CLE1-G1-50-18-180	500469	For gastroscope, 1 brush	5,0 mm	180 cm	50
⊗ CLE1-G2-50-18-180	500471	For gastroscope, 2 brushes	5,0 mm	180 cm	50
⊗ CLE1-C1-50-18-230	500460	For colonoscope, 1 brush	5,0 mm	230 cm	50
⊗ CLE1-C1-65-18-230	500462	For colonoscope, 1 brush	6,5 mm	230 cm	50
⊗ CLE1-C2-50-18-230	500464	For colonoscope, 2 brushes	5,0 mm	230 cm	50
⊗ CLE1-C2-30-16-280	500463	For colonoscope, 2 brushes	3,0 mm	280 cm	10
⊗ CLE1-V2-016	500480	For endoscope valve intake, 2 brushes	11,0 mm 5,0 mm	16 cm	10
⊗ CLE1-V2-016-50	502339	For endoscope valve intake, 2 brushes	11,0 mm 5,0 mm	16 cm	50
⊗ CLE1-C4-18-230	502496	For colonoscope, 1 brush, distal protection ball	5,0 mm	230 cm	50
⊗ CLE1-C6-18-230	502497	For colonoscope, 2 brushes, distal protection ball	5,0 mm	230 cm	50
⊗ CLE1-D2-18-230	502498	For colonoscope and endoscope valve intake, 2 brushes, distal protection ball	10,0 mm 5,0 mm	230 cm	50